



A High-frequency Ejector Delivers Top Performance

The Babylog VN500 ventilator was specially designed for premature babies and offers powerful high-frequency oscillation in addition to conventional ventilation. An important component of the ventilator is the expiratory valve, through which the patient exhales. The Infinity ID antenna module **1** transmits data from an Infinity ID hose system through the angled connector **2** to the ventilator, where the data are processed. Up to 18 liters of compressed gas at a pressure of a maximum of 2 bars is propelled into the ejector **3** each minute through two channels **4** measuring only 0.65 millimeters in width. The ejector itself is made of nickel silver, where the gas is fed into connector **5**. A lower pres-

sure, resulting from the ejection supports the “active expiration” of the patient, which is controlled via the silicone membrane **6**. The silicone membrane has a soft flat surface enabling a leak-free seal with minimal back pressure. A monocrystalline nickel disc vulcanized onto the back side ensures that it is perfectly flat.

The check membrane **7** prevents pendular breathing in the event of a possible failure of the device. The muffer **8** uses turbulence to reduce noise. The potential condensation is collected in the water trap **9**. The entire assembly can be cleaned manually or processed in an autoclave at 134 degrees Celsius.