

## The TransWarmer® Mattress

Critical warmth for the tiniest patient –  
Only nature builds a better nest.



The TransWarmer Mattress is absolutely essential when it comes to protecting newborns from dangerous cold stress. In emergency situations, infants are quite vulnerable and at risk of hypothermia. To keep the cold and stress at bay, the TransWarmer Mattress provides up to two hours of critical warmth.

### EASY-TO-USE AND RELIABLE

The TransWarmer Mattress is both flexible and easy to use: just firmly press the built-in trigger to activate. Its patented WarmGel® takes a mere 90 seconds to reach a temperature of 38°C (100.4 °F). WarmGel is a safe and proven latent heat store and source of steady warmth. It needs no external power supply, ensuring independence and flexibility no matter what the situation.

### COMFORTABLE

In addition to its special thermal properties, WarmGel is soft and pliable, creating a perfect “nest” for the newborn. The mattress also absorbs vibrations, thus lowering the risk of cerebral haemorrhage. This in turn allows caregivers to fully focus on the needs of the patient while meeting the challenges of transport.

### SAFE

Patient safety during such times is always paramount. The gel is non-toxic and does not irritate the skin or eyes, even in direct contact. Because the TransWarmer is designed for single use, it minimizes the risk of cross-contamination. It is also radio-lucent and MRI-compatible.

Warning: do not use in combination with warming therapy devices. Risk of hyperthermia!



MT-3302-2005

The TransWarmer® Mattress

**TECHNICAL DATA**

Size	25 x 41 cm / 10" x 16"
Content	Sodium acetate, water, and thickening agent (all contents are non-toxic)
Operating temperature	38°C at 22°C ambient temperature
Reaches operating temperature	within 90 seconds at ambient temperatures between 18°C and 28°C
Use in transport unit, delivery room, nursery, NICU, OR and ER. (Do not use in combination with warming devices.)	

**ORDER INFORMATION**

TransWarmer Mattress (6 single-use mattresses per box)	MP01414
--	---------

**HEADQUARTERS**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

**REGION EUROPE CENTRAL  
AND EUROPE NORTH**

Dräger Medical GmbH  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION EUROPE SOUTH**

Dräger Médical S.A.S.  
Parc de Haute  
Technologie d'Antony 2  
25, rue Georges Besse  
92182 Antony Cedex, France  
Tel +33 1 46 11 56 00  
Fax +33 1 40 96 97 20  
[d1mfr-contact@draeger.com](mailto:d1mfr-contact@draeger.com)

**REGION MIDDLE EAST, AFRICA,  
CENTRAL AND SOUTH AMERICA**

Dräger Medical GmbH  
Branch Office Dubai  
Dubai Healthcare City, P.O. Box 505108  
Dubai, United Arab Emirates  
Tel + 971 436 24 762  
Fax + 971 436 24 761  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

**REGION ASIA / PACIFIC**

Dräger Medical  
South East Asia Pte Ltd  
25 International Business Park  
#04-27/29 German Centre  
Singapore 609916, Singapore  
Tel +65 6572 4388  
Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**REGION NORTH AMERICA**

Dräger Medical, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 215 721 5400  
Toll-free +1 800 437 2437  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

**Manufacturer:**

Dräger Medical GmbH  
23542 Lübeck, Germany  
The quality management system at  
Dräger Medical GmbH is certified  
according to ISO 13485, ISO 9001  
and Annex II.3 of Directive 93/42/EEC  
(Medical devices).