

Dräger Blood Pressure Cuffs for Adult Patients

Our completely new line of Dräger blood pressure cuffs was designed for increased patient comfort, high accuracy and improved ease of use in any clinical setting.



NEW ERGONOMIC DESIGN

The new line of Dräger single patient use and reusable blood pressure cuffs offers an improved design for increased patient comfort and ease of use. The fastening system was redesigned to prevent skin irritation and provide improved flexibility.

Now available in an extended range of colour-coded sizes, they provide an excellent fit – regardless of patient size or age.

A COMPLETE PORTFOLIO

Dräger cuffs come in 9 sizes – from infant up to large adult – and include a specially marked extended range; a support for obese patients.


Regardless of the application, you can find a Dräger blood pressure cuff that can suit your patient needs. These quality cuffs can help you increase your efficiency, support hospital hygiene and improve patient comfort.



Single patient use cuffs



Reusable cuffs



ENHANCED PATIENT COMFORT

Our new cuff design features a more anatomical approach for both upper arm and thigh placement – with soft, rounded edges for increased patient comfort and easier handling for caregivers.

The hook-and-loop fastening system is positioned to avoid skin contact and irritation - making the cuffs more flexible and comfortable.

IMPROVED PATIENT SAFETY

Accurate blood pressure measurement depends on correct cuff fit. Dräger offers a choice of different cuff sizes for your patients. Dräger cuffs use color-coding and clearly visible, easy-to-understand symbols to support the correct cuff choice. Dräger cuffs are tested and validated for use with our Dräger Infinity® Patient Monitoring Systems.

Artery, index and range markings are visible on both sides of the cuff to ensure correct positioning. A newly developed Dräger connector, which can be easily opened using one hand, effectively prevents connection errors.

Single patient use cuffs Reusable cuffs

FOR HYGIENIC APPLICATIONS

Our single patient use cuffs support hospital hygiene and workflow by providing a disposable cuff solution for applications where the risk of infection or contamination is a factor. These cuffs are ideal for use in emergency wards, operating rooms, and infection or burn units. Because the hook-and-loop fastener attaches directly to the velour surface, the cuffs can be adjusted for maximum comfort and a best possible fit.

- Made of durable polyurethane
- Soft velour surface
- Free of Latex, silicone and PVC
- Soft, curved edges

FOR CONVENTIONAL APPLICATIONS

For conventional, everyday use – such as in regular wards – we recommend our line of reusable cuffs. Durable and easy to clean, these cuffs are designed to withstand the rigors of everyday use. Covered with soft, moisture-repellent material, the cuffs feel pleasant to the skin.

- Made of durable polyurethane
- Free of Latex, silicone and PVC
- Soft, curved edges
- Washable and submersible for easy disinfecting

The new NiBP cuff portfolio is compatible with the current Dräger Infinity® Patient Monitoring Systems.

D-253326-2009



Dräger connection system for adult cuffs

Our easy-to-use connector system supports one-handed use and is designed to prevent connection errors. To simplify handling, the connectors for neonatal and adult cuffs have different dimensions.

D-253336-2009



Free of Latex, silicone and PVC
The complete Dräger cuff portfolio is latex-, silicone- and PVC-free.

D-253335-2009



Dräger sizing information

Big cuff size markings are readable from a distance.

D-253411-2009



Dräger sizing information

To help you choose the correct cuff size, a symbol with the relevant patient dimensions is given on the patient side of every cuff.

ORDER INFORMATION

Single-patient use cuffs			Size	Part number	
NBP Cuff XXS	disp.	box of 10	Infant	8-13 / 13 cm	MP00924
NBP Cuff XS	disp.	box of 10	Child	12-19 / 19 cm	MP00925
NBP Cuff S	disp.	box of 10	Small Adult	17-25 / 29 cm	MP00926
NBP Cuff M	disp.	box of 10	Adult	23-33 / 33 cm	MP00928
NBP Cuff M+	disp.	box of 10	Adult Long	23-33 / 43 cm	MP00929
NBP Cuff M++	disp.	box of 10	Adult Extra	23-33 / 53 cm	MP00930
NBP Cuff L	disp.	box of 10	Large Adult	31-40 / 40 cm	MP00931
NBP Cuff L+	disp.	box of 10	Large Adult Long	31-40 / 55 cm	MP00932
NBP Cuff XL	disp.	box of 10	Thigh	38-50 / 50 cm	MP00934

Reusable cuffs			Size	Part number	
NBP Cuff XXS	reuse		Infant	8-13 / 13 cm	MP00911
NBP Cuff XS	reuse		Child	12-19 / 19 cm	MP00912
NBP Cuff S	reuse		Small Adult	17-25 / 29 cm	MP00913
NBP Cuff M	reuse		Adult	23-33 / 33 cm	MP00915
NBP Cuff M+	reuse		Adult Long	23-33 / 43 cm	MP00916
NBP Cuff M++	reuse		Adult Extra	23-33 / 53 cm	MP00917
NBP Cuff L	reuse		Large Adult	31-40 / 40 cm	MP00918
NBP Cuff L+	reuse		Large Adult Long	31-40 / 55 cm	MP00919
NBP Cuff XL	reuse		Thigh	38-50 / 50 cm	MP00921

Connector Hoses and Adapter		Size	Part number
NBP Connecting hose, adult		3,7 m	MP00953
Adapter old hose (1275275) - new connector, box of 10			MP00945

HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Effective 2010-09:

Dräger Medical AG & Co. KG
changes to
Dräger Medical GmbH

REGION EUROPE CENTRAL
AND EUROPE NORTH

Dräger Medical AG & Co. KG
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA,
CENTRAL AND SOUTH AMERICA

Dräger Medical AG & Co. KG
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel + 971 436 24 762
Fax + 971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

Manufacturer:

Dräger Medical AG & Co. KG
23542 Lübeck, Germany
The quality management system at
Dräger Medical AG & Co. KG is
certified according to ISO 13485,
ISO 9001 and Annex II.3 of Directive
93/42/EEC (Medical devices).